


HASA OPTIX
Belgium

General Catalogue 2024

Our innovative pre-packed pouches relevant to each ophthalmic procedure, help increase efficiency, convenience and peace of mind.



Availability of products

This catalogue may contain information on worldwide products and services, not all of which are available in every country. A reference to any product in this catalogue does not imply that such product or service is or will be available in your location. Any product or service referred to in this catalogue may be subject to different regulatory requirements depending on the country of use. Contact your local HASA OPTIX representative to learn more about the our products and services in your region.

General information and IP protection

© HASA OPTIX 2023. All rights reserved.

This catalogue and its contents are protected by copyright and intellectual property laws. Reproduction, distribution, or any unauthorized use of the material herein is strictly prohibited and may result in legal action.

A genuine desire to empower all those involved in ophthalmic surgery

HASA OPTIX holds a genuine desire to empower all those involved in ophthalmic surgery - with precision performance, procurement reliability and patient quality of care always at the forefront of our approach.

We design and manufacture premium, single-use, and fully recyclable ophthalmic instruments for both anterior and posterior segment surgeries.

Made of stainless steel, our exceptional precision instruments are supplied in recyclable pre-packed pouches relevant to each ophthalmic procedure. This not only ensures all surgeons have peace of mind regards sterilization and reliability, but also provides a cost controlled, simplified procurement and repository solution.

We are passionate about eliminating cross contamination and reducing risk management whilst increasing efficiency, convenience, and reliability.

Our innovative solution provides total peace of mind, confidence, and satisfaction.



What Sets Us Apart?

Quality improvement and cost-savings are of great importance in healthcare in general and in ophthalmic surgery. Hospitals, more than ever, must face the challenges of working profitably and obtaining optimum quality at the same time. In a changing world, we believe there are new ways that will allow to facilitate work of all stakeholders involved in ophthalmic surgery, allowing a higher pharmacoeconomic value.



PREMIUM QUALITY

Premium quality disposable recyclable instruments offering the same constant quality as reusable instruments but with increased convenience and efficiency.



SIMPLIFY

Our pre-packed pouches help reduce the number of SKUs, contributing towards a simpler, more efficient repository solution.



COST CONTROL

HASA OPTIX instruments provide the ideal procurement solution by aiding the precise calculation of surgery costs.



PEACE OF MIND

Our innovative solution provides total peace of mind ensuring the correct surgical instrumentals are easily accessible.



INFECTION CONTROL

Our single use and pre-sterilized instruments ensure all surgeons can be confident with their patient safety of care.



ENVIRONMENTALLY FRIENDLY

Supplied in recyclable pre-packed pouches relevant to each ophthalmic procedure, helping to reduce waste and environmental impact.

Important Information:

Individual instruments are supplied in a box of 10 pcs, unless otherwise specified.

Instrument sets are supplied in a box of 10 sets, unless otherwise specified.

When **MOQ** is indicated, it signifies that the MOQ for that SKU is 1000 pieces.

FORCEPS

NEEDLE HOLDERS

MANIPULATORS, CHOPPERS & SPATULAS

PROBES & DILATORS

MARKERS & GAUGES

SCISSORS


MICRO INSTRUMENTS

SPECULUMS

CURETTES & RETRACTORS

DEPRESSORS & HOOKS

MISCELLANEOUS



FORCEPS

FORCEPS

Capsulorhexis Forceps

All Capsulorhexis Forceps are available for incision from 1.8mm to 3.2mm

Capsulorhexis Forceps, Curved Thin Shafts,
Flat Cross-Action Handle

651000



Capsulorhexis Forceps, 45° Angled Thin Shafts,
Flat Cross-Action Handle

651001



Capsulorhexis Forceps, Curved Thin Shafts, Cross-Action,
Round Knurled Handle

651002




Capsulorhexis Forceps, 45° Angled Thin Shafts, Cross-Action,
Round Knurled Handle

651003



Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle

651141



FORCEPS

Capsulorhexis Forceps

All Capsulorhexis Forceps are available for incision from 1.8mm to 3.2mm


Capsulorhexis Forceps, Curved Thin Shafts,
Cross-Action Ripped Handle

651143



Capsulorhexis Forceps, 45° Angled Thin Shafts,
Round Knurled Handle

651147




Capsulorhexis Forceps, 45° Angled Thin Shafts,
Cross-Action Ripped Handle

651149




Capsulorhexis Forceps, 45° Angled Thin Shafts,
Flat Handle

651150



Capsulorhexis Forceps, Curved Thin Shafts,
Flat Handle

651153



FORCEPS

Fixation Forceps

Toothed Forceps, Straight Shafts With 5mm Tying Platforms, Flat Handle, 1x2 Teeth

551000

| SKU | Size |
|--------|---------|
| 551000 | 0.10 mm |
| 551001 | 0.12 mm |
| 551002 | 0.15 mm |
| 551003 | 0.30 mm |
| 551004 | 0.40 mm |
| 551005 | 0.50 mm |



Toothed Forceps, Straight Shafts With 5mm Tying Platforms, Flat Handle, 1x2 Teeth

551006

| SKU | Size |
|--------|---------|
| 551006 | 0.10 mm |
| 551007 | 0.12 mm |
| 551008 | 0.15 mm |
| 551009 | 0.30 mm |
| 551010 | 0.40 mm |
| 551011 | 0.50 mm |



Toothed Forceps, Straight Shafts With 5mm Tying Platforms, Flat Handle, 1x2 Teeth

551012


| SKU | Size |
|--------|---------|
| 551012 | 0.10 mm |
| 551013 | 0.12 mm |
| 551014 | 0.15 mm |
| 551015 | 0.30 mm |
| 551016 | 0.40 mm |
| 551017 | 0.50 mm |



Toothed Forceps, Straight Shafts With 5mm Tying Platforms, Round Knurled Handle, 1x2 Teeth

551018

| SKU | Size |
|--------|---------|
| 551018 | 0.10 mm |
| 551019 | 0.12 mm |
| 551020 | 0.15 mm |
| 551021 | 0.30 mm |
| 551022 | 0.40 mm |
| 551023 | 0.50 mm |



FORCEPS

Fixation Forceps

Toothed Forceps, Straight Shafts With 5mm Tying Platforms, Castroviejo Handle, 1x2 Teeth

551024

| SKU | Size |
|--------|---------|
| 551024 | 0.10 mm |
| 551025 | 0.12 mm |
| 551026 | 0.15 mm |
| 551027 | 0.30 mm |
| 551028 | 0.40 mm |
| 551029 | 0.50 mm |



Curved Toothed Forceps, Curved Shafts With 5mm Tying Platforms, Flat Handle, 1x2 Teeth

551037

| SKU | Size |
|--------|---------|
| 551037 | 0.10 mm |
| 551038 | 0.12 mm |
| 551039 | 0.15 mm |
| 551040 | 0.30 mm |
| 551041 | 0.40 mm |
| 551042 | 0.50 mm |



Colibri Conjunctiva Toothed Forceps, Curved Shafts With 5mm Tying Platforms, Flat Handle, 1x2 Teeth

551043

| SKU | Size |
|--------|---------|
| 551043 | 0.10 mm |
| 551044 | 0.12 mm |
| 551045 | 0.15 mm |
| 551046 | 0.30 mm |
| 551047 | 0.40 mm |
| 551048 | 0.50 mm |



Colibri Toothed Forceps, Curved Shafts With 5mm Tying Platforms, Flat Handle, 1x2 Teeth

551049

| SKU | Size |
|--------|---------|
| 551049 | 0.10 mm |
| 551050 | 0.12 mm |
| 551051 | 0.15 mm |
| 551052 | 0.30 mm |
| 551053 | 0.40 mm |
| 551054 | 0.50 mm |



FORCEPS

Fixation Forceps

Colibri Toothed Forceps, Curved Shafts
With 5mm Tying Platforms, Flat Handle, 1×2 Teeth

551055


| SKU | Size |
|--------|---------|
| 551055 | 0.10 mm |
| 551056 | 0.12 mm |
| 551057 | 0.15 mm |
| 551058 | 0.30 mm |
| 551059 | 0.40 mm |
| 551060 | 0.50 mm |



Curved Toothed Forceps, Curved Shafts
With 5mm Tying Platforms, Round Knurled Handle, 1×2 Teeth

551061

| SKU | Size |
|--------|---------|
| 551061 | 0.10 mm |
| 551062 | 0.12 mm |
| 551063 | 0.15 mm |
| 551064 | 0.30 mm |
| 551065 | 0.40 mm |
| 551066 | 0.50 mm |



Curved Toothed Forceps, Curved Shafts
With 5mm Tying Platforms, Castroviejo Handle, 1×2 Teeth

551067

| SKU | Size |
|--------|---------|
| 551067 | 0.10 mm |
| 551068 | 0.12 mm |
| 551069 | 0.15 mm |
| 551070 | 0.30 mm |
| 551071 | 0.40 mm |
| 551072 | 0.50 mm |



Kelman-McPherson Forceps, Angled Shafts
With 5mm Tying Platforms, Flat Handle, 1×2 Teeth

551080

| SKU | Size |
|--------|---------|
| 551080 | 0.10 mm |
| 551081 | 0.12 mm |
| 551082 | 0.15 mm |
| 551083 | 0.30 mm |
| 551084 | 0.40 mm |
| 551085 | 0.50 mm |




FORCEPS

Fixation Forceps

Kelman-McPherson Forceps, Angled Shafts
With 5mm Tying Platforms, Flat Handle, 1×2 Teeth

551086

| SKU | Size |
|--------|---------|
| 551086 | 0.10 mm |
| 551087 | 0.12 mm |
| 551088 | 0.15 mm |
| 551089 | 0.30 mm |
| 551090 | 0.40 mm |
| 551091 | 0.50 mm |



Kelman-McPherson Forceps, Angled Shafts
With 5mm Tying Platforms, Flat Handle, 1×2 Teeth

551092

| SKU | Size |
|--------|---------|
| 551092 | 0.10 mm |
| 551093 | 0.12 mm |
| 551094 | 0.15 mm |
| 551095 | 0.30 mm |
| 551096 | 0.40 mm |
| 551097 | 0.50 mm |



Kelman-McPherson Forceps, Angled Shafts
With 5mm Tying Platforms, Round Knurled Handle, 1×2 Teeth

551098

| SKU | Size |
|--------|---------|
| 551098 | 0.10 mm |
| 551099 | 0.12 mm |
| 551100 | 0.15 mm |
| 551101 | 0.30 mm |
| 551102 | 0.40 mm |
| 551103 | 0.50 mm |



Kelman-McPherson Forceps, Angled Shafts
With 5mm Tying Platforms, Castroviejo Handle, 1×2 Teeth

551104

| SKU | Size |
|--------|---------|
| 551104 | 0.10 mm |
| 551105 | 0.12 mm |
| 551106 | 0.15 mm |
| 551107 | 0.30 mm |
| 551108 | 0.40 mm |
| 551109 | 0.50 mm |



FORCEPS

Fixation Forceps

Polack Corneal Suturing Forceps, Curved Shafts
Without Platforms, Flat Handle, 0.12mm, 1×2 Teeth

551119



Notched Forceps, Straight Shafts With Tying Platforms,
Flat Handle, 0.25mm Notched


551122

| SKU | Overall length |
|--------|----------------|
| 551122 | 85 mm |
| 551123 | 105 mm |
| 551124 | 115 mm |



Notched Forceps, Straight Shafts With Tying Platforms,
Round Knurled Handle, 0.25mm Notched

551125



Notched Forceps, Straight Shafts With Tying Platforms,
Castroviejo Handle, 0.25mm Notched

551126



FORCEPS

Fixation Forceps

Notched Forceps, Curved Shafts With Tying Platforms, Flat Handle, 0.25mm Notched

551128

| SKU | Overall length |
|--------|----------------|
| 551128 | 85 mm |
| 551129 | 105 mm |
| 551130 | 115 mm |




Pierse-Colibri Corneal Forceps, Curved Shafts With Tying Platforms, Flat Serrated Handle, 0.25mm Notched

551131



Notched Forceps, Curved Shafts With Tying Platforms, Round Handle, 0.25mm Notched

551132



Notched Forceps, 45° Angled Straight Shafts With Tying Platform, Flat Handle, 0.25mm Notched

551135

| SKU | Overall length |
|--------|----------------|
| 551135 | 85 mm |
| 551136 | 105 mm |
| 551137 | 115 mm |




FORCEPS

Fixation Forceps

Notched Forceps, 45° Angled Straight Shafts
With Tying Platform, Round Handle, 0.25mm Notched

551138




Notched Forceps, 45° Angled Straight Shafts
With Tying Platform, Castroviejo Handle, 0.25mm Notched

551139



Hofmann-Polack Corneal Suturing Forceps,
45° Angled Straight Shafts Without Platform, Flat Handle,
0.25mm Notched

551141



Colibri Toothed Forceps, Curved Shafts, Serrated Handle,
0.12mm, 1x2 Teeth

551187




FORCEPS

Fixation Forceps

Paufigue Forceps, Oblique Teeth 0.7mm And Tying Platforms

651080



Corneal Forceps, Oblique Teeth, Straight With Tying Platforms, Dipped Serrated Handle

651081

| SKU | Size |
|--------|---------|
| 651081 | 0.12 mm |
| 651082 | 0.20 mm |



Corneal Forceps, Oblique Teeth 0.12mm, Curved With Tying Platforms, Dipped Serrated Handle

651085



FORCEPS

Tissue Forceps


Adson Forceps, Straight Shafts Without Platforms, Swallow-Tail Handle, 0.40mm, 1x2 Teeth

551120



Moorfields Forceps, Blunt Tipped Forceps

551235



Gillis Tissue Forceps

551240

| SKU | Teeth |
|--------|---------|
| 551240 | 0.50 mm |
| 551241 | 0.40 mm |
| 551242 | 0.30 mm |



Adson Forceps Non-Toothed, Serrated Jaws

651120



FORCEPS

Tying Forceps

Tying Forceps, Straight Shafts With 4.5mm Tying Platform, Flat Handle


551169

| SKU | Overall length |
|--------|----------------|
| 551169 | 85 mm |
| 551170 | 105 mm |
| 551171 | 115 mm |



Tying Forceps, Straight Shafts With 4.5mm Tying Platform, Round Knurled Handle

551172



Tying Forceps, Curved Shafts With 6mm Tying Platform, Flat Handle


551175

| SKU | Overall length |
|--------|----------------|
| 551175 | 85 mm |
| 551176 | 105 mm |
| 551177 | 115 mm |



Tying Forceps, Curved Shafts With 6mm Tying Platform, Round Knurled Handle

551178



FORCEPS

Tying Forceps

Tying Forceps, Curved Shafts With 6mm Tying Platform, Castroviejo Handle

551179



Tying Forceps, 45° Angled Shafts With 4.5mm Tying Platform, Flat Handle



551182

| SKU | Overall length |
|--------|----------------|
| 551181 | 85 mm |
| 551182 | 105 mm |
| 551183 | 115 mm |



Tying Forceps, 45° Angled Shafts With 4.5mm Tying Platform, Round Knurled Handle

551184



Tying Forceps, 45° Angled Shafts With 4.5mm Tying Platform, Castroviejo Handle

551185



FORCEPS

Tying Forceps

McPherson Tying Forceps, With 7.5mm Tying Platforms

651187



FORCEPS

Lens Forceps

Lens Forceps, Curved Shafts With 90° Angled 1.0mm Tips, Cross-Action Handle

551188



Lens Forceps, Curved Shafts With 1.0mm Ring Tips, Cross-Action Handle

551189



Lens Forceps, Angled Shafts With 1.0mm Round Tips, Flat Handle

551190



Lens Forceps, 45° Angled Shafts With Smooth Flat Jaws, Cross-Action Handle

551191



FORCEPS

Lens Forceps

Kelman-McPherson Lens Forceps, Angled Shafts With Smooth Jaws, Castroviejo Handle

551192

| SKU | Overall Length |
|--------|----------------|
| 551192 | 102mm |
| 551193 | 105mm |



Lens Implantation Forceps, Smooth And Curved Jaws With Lock, Cross-Action Handle


551194

| SKU | Lock |
|--------|---------|
| 551194 | with |
| 551195 | without |



Kelman-McPherson Lens Forceps, Angled Shafts, Flat Handle

551199



ICL Cartridge Loading Forceps, Cross-Action

651071



FORCEPS

Lens Forceps

ICL Cartridge Loading Forceps, Flat Handle

651073



IOL Forceps

651076



Dennis Forceps, Straight Shafts, Flat Handle

MOQ

651086



Lenticule Extraction Forceps

651160




FORCEPS

Lens Forceps

Lenticule Removal Forceps, Serrated Jaws,
Angled, Castroviejo Handle

651161



FORCEPS

LASIK Forceps


Perone LASIK Flap Forceps, 2mm Rings Tips

551211



Perone LASIK Flap Forceps, 2mm Long Tips

551212



LASIK Retreatment Forceps, 1x2Teeth, 0.50mm

551225



LASIK Multipurpose Forceps, Curved Jaws

551230



FORCEPS

Prechopper Forceps

Akahoshi Prechopper Forceps, Flat Handle, Straight Shafts And Cross-Action, Maximum Tip Opening 3mm

551200



Akahoshi Prechopper Forceps, Round Handle, Straight Shafts And Cross-Action

551201




Akahoshi Prechopper Forceps, Cross-Action, 11mm Angled Long

551205



Akahoshi Prechopper Forceps, Cross-Action, 10mm Angled Shafts

551206




FORCEPS

Prechopper Forceps


Akahoshi Prechopper Forceps, Cross-Action, 12mm Angled Shafts

551207



Akahoshi Prechopper Forceps, Cross-Action, 12mm Curved Shafts

551208



FORCEPS

Chalazion Forceps

Chalazion Forceps, Serrated Handle
With Locking Thumb Screw, Round 8mm

551213



Chalazion Forceps, Serrated Handle
With Locking Thumb Screw, Oval 8×10mm

551214



Chalazion Forceps, Serrated Handle
With Locking Thumb Screw, Round 15mm

551215



Chalazion Forceps, Serrated Handle
With Locking Thumb Screw, Round 12mm

551216



FORCEPS

Chalazion Forceps

Chalazion Forceps, Serrated Handle
With Locking Thumb Screw, Oval 12×14mm

551217



Ciliary Forceps, Flat Handle

551218



Chalazion Forceps, 11mm Diameter, Round Plate

551236



Chalazion Forceps, Oval Plate, 21×12mm

551237




FORCEPS

General Forceps

Scleral Plug Forceps, Swallow-Tail Handle, Cross-Action Shafts

551219



Jeweler Style Forceps, Straight, 0.05x0.01mm Tip

551220



Jeweler Style Forceps, Angled 45°, Extra Delicate, 0.05x0.01mm Tip

551222



Jeweler Style Forceps, Curved, 0.07x0.03mm Tip

551223



FORCEPS

General Forceps

Landers Vitrectomy Lens Forceps,
Shafts Angled 45° With Curved Grooved Blades

551224



Mosquito Forceps, Straight 24mm Serrated Jaws

552034



Mosquito Forceps, Curved 24mm Serrated Jaws

552035



Towel Clamps, Flat Handle

552042




FORCEPS

General Forceps

Meibomian Forceps


651070



Florakis Endothial Forceps

MOQ

651090



Iris Claw Forceps

651091



Mermoud Forceps

651156




FORCEPS

General Forceps

Plain Forceps

651162





NEEDLE HOLDERS



NEEDLE HOLDERS

Needle Holder, 9mm Straight Jaws, Flat Handle

552001

| SKU | Lock |
|--------|---------|
| 552000 | with |
| 552001 | without |



Needle Holder, 11mm Straight Jaws, Flat Handle

552002

| SKU | Lock |
|--------|---------|
| 552002 | with |
| 552003 | without |



Needle Holder, 6mm Straight Jaws, Round Knurled Handle

552004

| SKU | Lock |
|--------|---------|
| 552004 | with |
| 552005 | without |



Needle Holder, 9mm Straight Jaws, Flat Handle

552006

| SKU | Lock |
|--------|---------|
| 552006 | with |
| 552007 | without |



Needle Holder, 11mm Curved Jaws, Round Knurled Handle

552008

| SKU | Lock |
|--------|---------|
| 552008 | with |
| 552009 | without |



NEEDLE HOLDERS

Needle Holder, 9mm Curved Jaws, Round Knurled Handle

552010

| SKU | Lock |
|--------|---------|
| 552010 | with |
| 552011 | without |



Needle Holder, 9mm Curved Jaws, Flat Handle

552012

| SKU | Lock |
|--------|---------|
| 552012 | with |
| 552013 | without |



Needle Holder, 11mm Curved Jaws, Flat Handle

552014

| SKU | Lock |
|--------|---------|
| 552014 | with |
| 552015 | without |



Needle Holder, 9mm Curved Jaws, Flat Handle

552016


| SKU | Lock |
|--------|---------|
| 552016 | with |
| 552017 | without |



Needle Holder, 9mm Curved Jaws, Round Knurled Handle

552018

| SKU | Lock |
|--------|---------|
| 552018 | with |
| 552019 | without |



NEEDLE HOLDERS

Needle Holder, 9mm Curved Jaws, Round Knurled Handle

552020


| SKU | Lock |
|--------|---------|
| 552020 | with |
| 552021 | without |



Needle Holder, 10mm Curved Jaws, Round Knurled Handle

552022

| SKU | Lock |
|--------|---------|
| 552022 | with |
| 552023 | without |



Needle Holder, 9mm Curved Jaws, Round Knurled Handle

552024


| SKU | Lock |
|--------|---------|
| 552024 | with |
| 552025 | without |



Needle Holder, 6mm Curved Jaws, Round Knurled Handle

552026

| SKU | Lock |
|--------|---------|
| 552026 | with |
| 552027 | without |



Needle Holder, 9mm Curved Jaws, Flat Handle

552028

| SKU | Lock |
|--------|---------|
| 552028 | with |
| 552029 | without |



NEEDLE HOLDERS

Needle Holder, 9mm Curved Jaws, Round Knurled Handle

552030

| SKU | Lock |
|--------|---------|
| 552030 | without |



Needle Holder, 9mm Curved Jaws, Flat Handle

552032


| SKU | Lock |
|--------|---------|
| 552032 | without |



Needle Holder, 9mm Curved Jaws, Round Knurled Handle

652021

| SKU | Lock |
|--------|---------|
| 652021 | without |





MANIPULATORS CHOPPERS & SPATULAS



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

Kuglen Iris Hook, K-Shaped Tip

554000

| SKU | Tip |
|--------|----------|
| 554000 | Angled |
| 554001 | Straight |



Jaffe-Knolle Iris Hook

554002

| SKU | Tip |
|--------|----------|
| 554002 | Angled |
| 554003 | Straight |



Pupil Snapper Hook,
With Smoothly Polished Notch & Sharp Inner Edge

554004

| SKU | Tip |
|--------|--------|
| 554004 | Angled |



Lens Manipulator, T-Shaped Tip

554006

| SKU | Tip |
|--------|----------|
| 554005 | Straight |
| 554006 | Angled |



Lens Manipulator, Y Nucleus Rotator

554007

| SKU | Tip |
|--------|------------|
| 554007 | Vertical |
| 554087 | Horizontal |



Lens Manipulator, V-Shaped Tip

554008

| SKU | Tip |
|--------|----------|
| 554008 | Angled |
| 554009 | Straight |



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

Lens Manipulator, 0.7mm Tip Length, Angled

554010



Lens Manipulator, 1.1mm Tip Length, Curved

554011



Nucleus Manipulators, Paddle Shaped Tip

554012



Nucleus Mushroom Manipulators

554013



Nucleus Manipulators, Connor Wand, Straight

554014



Nucleus Manipulators, Connor Wand, Angled, Round

554015



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

Nucleus Manipulators, Connor Wand, Angled, Oval

554016



Iris Hook, H-Shaped Tip, Angled

554018



Repositor, Semi-Sharp Tip, 0.7mm

554027



Repositor, Semi-Sharp Tip, 0.5mm

554028



Kansas Nucleus Bisector

554053



Phaco Chopper

554054

| SKU | Hand |
|--------|-------|
| 554054 | Left |
| 554055 | Right |



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

Phaco Chopper, "M" Hook

554056



Phaco Chopper, Ambidextrous Nucleus Chopper

554057



Phaco Chopper, Angled Nucleus Chopper

554058



Phaco Chopper, Curved Chopper With Two Faceted Chopping Surfaces

554059



Phaco Chopper, 1.8mm Tip

554061

| SKU | Blade |
|--------|----------|
| 554060 | Straight |
| 554061 | Curved |



Phaco Chopper, Angle Shaft, 1.6mm Tip

554062

| SKU | Hand |
|--------|-------|
| 554062 | Left |
| 554063 | Right |



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

Haefliger Phaco Chopper (Chopper & Manipulator in One)

554064



Phaco Chopper, Special Inner Blade Is Offset 45°

554065



Kuglen Iris Hook, "S" Shaped 0.6mm Hook

554066



Kuglen Iris Hook, Curved Hook

554067




Intra-LASIK Flap Separator, Round, Sharp Tipped Separator

554071



Corneal Epithelium Spatula

554082



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

Epi-Peeler

554083



Hoe Type Epi-Peeler, Square Shaped

554084



Koch Stop & Chop Manipulator

554085



Rosen Phaco Chopper

554088



Koch Phaco Spatula

554095



LASIK Flap Lifter/Elevator Binder

554504



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

ICL Manipulator

654073

| SKU | Shaft | Tip diameter |
|--------|--------|--------------|
| 654073 | Angled | 1.0 mm |
| 654074 | Curved | 1.0 mm |
| 654075 | Curved | 0.6 mm |
| 654076 | Curved | 0.6 mm |



Lund Iris Hook, Angled 45°

654088



Manipulator, Lester Angled

654097



Sinskey Hook, Straight

654098



Lebuisson Manipulator

654107



Manipulator Cylinder

654112



MANIPULATORS, CHOPPERS & SPATULAS

Single-Ended

Manipulator Govan

654113



Denis Nucleus Chopper

654114



Manipulator Clayman Spatula, Flat

654119



Manipulator Sinskey, Tip 0.2mm

654120

| SKU | Altitude |
|--------|----------|
| 654120 | 0.7 mm |
| 654121 | 1.1 mm |



Multifunctional Manipulator

654128



Paufique Spatula

654129




MANIPULATORS, CHOPPERS & SPATULAS

Double-Ended


Repositor, Double Ended Iris Spatula

554029




Double Ended Drysdale & Chopper

554501



Fukasaku & Sinsky

554502




MANIPULATORS, CHOPPERS & SPATULAS

Double-Ended


LASIK Flap Separator/Elevator

554503




Double-Ended Manipulator,
Drysdale/H-Shaped Tip

654085



Double-Ended Manipulator,
Drysdale/Y Pusher

654086




MANIPULATORS, CHOPPERS & SPATULAS

Double-Ended


Double-Ended Manipulator,
Cylinder/Y Pusher

654087




Double-Ended Manipulator,
Drysdale/ Fukasaku Small Pupil Snapper

654089



Double-Ended Manipulator,
Agarwal Phaco Chopper/Fukasaku Small Pupil Snapper

654090




MANIPULATORS, CHOPPERS & SPATULAS

Double-Ended


Double-Ended Manipulator,
Drysdale/Koch Chopper

654096




Double-Ended Manipulator,
Drysdale/Curved Chopper

654100



Double-Ended Manipulator,
Fukasaku Small Pupil Snapper/Cylinder

654104




MANIPULATORS, CHOPPERS & SPATULAS

Double-Ended

Double-Ended Manipulator,
Koch Chop/Chopper

654106




Lenticule Separator/Elevator

654108



Double-Ended Manipulator,
Sinskey/Lester

654109




MANIPULATORS, CHOPPERS & SPATULAS

Double-Ended


Double-Ended Manipulator,
Lund Iris Hook/Phaco Chopper

654110




Double-Ended Manipulator,
Mushroom/Chopper

654111



Double-Ended Manipulator,
Sweep/Y

654118



MANIPULATORS, CHOPPERS & SPATULAS

Double-Ended

Double-Ended Reversed Sinskey/Sinskey

654123



LASIK Flap Lifter/Sinskey

654124





PROBES & DILATORS

PROBES & DILATORS

Lacrimal Probe, Stainless Steel Set of 4 pieces

731

| SKU | Size |
|-----|-----------|
| 731 | 0.5/0.6mm |
| | 0.7/0.8mm |
| | 0.9/1mm |
| | 1.1/1.2mm |



Lacrimal Probe, Pigtail Curved Probes

732

| SKU | Diameter |
|-----|----------|
| 732 | 6mm |
| 733 | 7mm |
| 734 | 8mm |
| 735 | 9mm |



Punctal Dilator

736

| SKU | Taper | Diameter |
|-----|-------|----------|
| 736 | 16mm | 0.3mm |
| 737 | 15mm | 0.4mm |
| 738 | 15mm | 0.3mm |



Punctal Dilator, Sharp 17mm Taper & Sharp 13mm Taper

739



Double Ended Dilator 15mm & 19mm

741



Lens Loop, 5x7mm

554034

| SKU | Loop |
|--------|-------|
| 554034 | 5x7mm |
| 554035 | 3x8mm |





MARKERS, CALIPERS & GAUGES

MARKERS & GAUGES

LASEK Alcohol Well

554036

| SKU | Diameter |
|--------|----------|
| 554036 | 7.0mm |
| 554037 | 8.0mm |
| 554038 | 9.0mm |
| 654039 | 9.5mm |



Optical Zone, Target Centre, 4.0mm Diameter

558000



CK LASIK Marker, Target Centre, 9.5mm Diameter

558001



Mendez Degree Gauge, Calibrated Every 10° From 0° to 180°, External Diameter 19.0mm, IØ 12.0mm

558002



Scleral Marker, Double Ended, Aphakic 3.5/4.0mm


558003



MARKERS & GAUGES

LASIK Corneal Marker,
Two Blunt Fins & Separated at 45°Angle

558004



Braunstein Fixed Caliper
Fine Tips Mark Sclerotomy Sites For Vitreous Surgery, Aphakic 3.5/4.0mm

558005



Braunstein Fixed Caliper,
Fine Tips Mark Sclerotomy Sites For Vitreous Surgery, Aphakic 3.5/4.0mm

558006



Castroviejo Caliper,
Scale Reads On One Side, Measure From 0mm To 20mm

558009





SCISSORS

SCISSORS

Castroviejo Scissors, Curved, Blunt Tips,
Tip To Pivot Length 14.0mm

553000



Castroviejo Scissors, Straight, Blunt Tips,
Tip To Pivot Length 14.0mm

553001



Conjunctiva Scissors, Curved, Blunt Tips

553002

| SKU | Tip to pivot length |
|--------|---------------------|
| 553002 | 11 mm |
| 553003 | 12 mm |
| 553004 | 16.5 mm |



Castroviejo Corneal Scissors, Curved, Blunt Tips,
Tip To Pivot Length 16mm

553006

| SKU | Blades |
|--------|--------------|
| 553005 | Large blades |
| 553006 | Small blades |



Castroviejo Corneal Scissors, Blunt Tips And Lower Blade,
Tip To Pivot Length 13mm

553007

| SKU | Blades |
|--------|--------|
| 553007 | Left |
| 553008 | Right |



SCISSORS

Castroviejo Corneal Scissors, Left, Blunt Tips And Lower Blade, Tip To Pivot Length 11mm

553009


| SKU | Blades |
|--------|--------|
| 553009 | Left |
| 553010 | Right |



Castroviejo Corneal Scissors, Curved, Blunt Tips

553011

| SKU | Tip to pivot length |
|--------|---------------------|
| 553011 | 12mm |
| 553012 | 16mm |



Castroviejo Corneal Scissors, Angled, Blunt Tips & Lower Blade, Tip To Pivot Length 13mm

553013



Vannas Scissors, Curved, Sharp Pointed Tips

553014

| SKU | Tip to pivot length |
|--------|---------------------|
| 553014 | 10mm |
| 553015 | 16mm |



Vannas Scissors, Straight, Sharp Pointed Tips Tip To Pivot Length 10mm

553016



SCISSORS

Vannas Scissors, Angled, Sharp Pointed Tips, Tip To Pivot Length 9.0mm

553017



Iris Scissors, Curved, Sharp Tips, Cutting Length 6.0mm, Blade Length 11mm

553018



Iris Scissors, Curved, Sharp Tips, Blade Length 13mm

553019



Iris Scissors, Straight, Sharp Tips, Cutting Length 6.0mm, Blade Length 11mm

553020



Iris Scissors, Straight, Sharp Tips, Blade Length 12mm

553021



SCISSORS

Westcott Stitch Scissors, Curved, Sharp Pointed Tips

553022

| SKU | Tip to pivot length |
|--------|---------------------|
| 553022 | 14mm |
| 553023 | 16mm |
| 553024 | 11mm |



Vannas Capsulotomy Scissors, Curved, Sharp Pointed Tips

553025

| SKU | Blade length |
|--------|--------------|
| 553025 | 11.5mm |
| 553026 | 15mm |



Vannas Capsulotomy Scissors, Straight, Sharp Pointed Tips

553027

| SKU | Blade length |
|--------|--------------|
| 553027 | 13mm |
| 553028 | 16mm |



Vannas Capsulotomy Scissors, Angled, Sharp Pointed Tips

553029

| SKU | Blade length |
|--------|--------------|
| 553029 | 12mm |
| 553030 | 14mm |



Tenotomy Scissors, Curved, sharp Pointed Tips

553031

| SKU | Tip to pivot length | Overall length |
|--------|---------------------|----------------|
| 553031 | 16mm | 115mm |
| 553032 | 20mm | 125mm |



SCISSORS

Westcott Stitch Scissors, Curved, Sharp Pointed Tips,
Tip To Pivot Length 16mm

553034



Vannas Capsulotomy Scissors, Angled, Shaft Pointed Tips,
Tip To Pivot Length 16mm

553035



Vannas Scissors, Straight, Sharp Pointed Tips,
Tip To Pivot Length 10mm

553036



Eye Scissors, Straight, Sharp Pointed

553038



Eye Scissors, Straight, Blunt Ended

553039




| SKU | Overall length |
|--------|----------------|
| 553039 | 90mm |
| 553040 | 130mm |

SCISSORS

Westcott Stitch Scissors, Curved, Sharp Tips, Tip To Pivot Length 21mm

653034



Iris Scissors Hasa Optix, Straight, Sharp

653041



Iris Scissors Hasa Optix, Curved, Sharp

653042



Iris Scissors Hasa Optix, Curved, Blunt

653043



Iris Scissors Hasa Optix, Straight, Blunt

653044



SCISSORS


Westcott Scissors, Straight, Blunt

653045



Westcott Stitch Scissors, Curved, Blunt, Tip To Pivot Length 22mm

653055




Westcott Stitch Scissors, Straight, Sharp, Tip To Pivot Length 24mm

653056




IOL Cutter

653058



Noyes Iris Scissors Straight Blunt

653060



SCISSORS

Noyes Iris Scissors Straight Sharp

653061




Noyes Iris Scissors Curved Blunt

653062



Noyes Iris Scissors Curved Sharp

653063



Stevens, Tenotomy Scissors, Short Blades, Blunt Tips, 34mm From Tip to Pivot, Straight

753013

MOQ



Enucleation Scissors, Curved Blades, 40mm Tips to Pivot, Strong Curve

753035

MOQ





MICRO INSTRUMENTS

MICRO INSTRUMENTS

Micro Forceps



Vitreoretinal Forceps,
Short End-Gripping Forceps (ILM)

556026

| SKU | Size |
|--------|------|
| 556026 | 23G |
| 656026 | 25G |
| 856026 | 27G |



Vitreoretinal Forceps, Gripping Forceps,
1.6mm Smooth, Rounded Tips

556028

| SKU | Size |
|--------|------|
| 556028 | 23G |
| 656028 | 25G |
| 856028 | 27G |



Vitreoretinal Forceps, Serrated, Straight

556030


| SKU | Size |
|--------|------|
| 556030 | 23G |
| 656030 | 25G |
| 856030 | 27G |



Vitreoretinal Forceps, Straight, With 1x2 Teeth

556032

| SKU | Size |
|--------|------|
| 556032 | 23G |
| 656032 | 25G |
| 856032 | 27G |



Vitreoretinal Forceps, Membrane Peeler, Smooth, 30°

556035

| SKU | Size |
|--------|------|
| 556035 | 23G |
| 656035 | 25G |
| 856035 | 27G |



MICRO INSTRUMENTS

Micro Forceps

Vitreoretinal Forceps, Membrane Peeler, Grooved, 30°

556037

| SKU | Size |
|--------|------|
| 556037 | 23G |
| 656037 | 25G |
| 856037 | 27G |



Vitreoretinal Forceps, Long End-Gripping Forceps

556038

| SKU | Size |
|--------|------|
| 556038 | 23G |
| 656038 | 25G |
| 856038 | 27G |



Vitreoretinal Forceps, End-Gripping Forceps Round

556039


| SKU | Size |
|--------|------|
| 556039 | 23G |
| 656039 | 25G |
| 856039 | 27G |



Vitreoretinal Forceps, Asymmetric End-Gripping

556042

| SKU | Size |
|--------|------|
| 556042 | 25G |
| 856042 | 27G |



Vitreoretinal Forceps, Pointed Short Forceps Flat

556043

| SKU | Size |
|--------|------|
| 556043 | 23G |
| 654043 | 25G |
| 856043 | 27G |



MICRO INSTRUMENTS

Micro Forceps

Vitreoretinal Forceps, Smooth Serrated

556044

| SKU | Size |
|--------|------|
| 556044 | 23G |
| 656044 | 25G |
| 856044 | 27G |



Vitreoretinal Wide End-gripping Forceps (ILM)

556050

| SKU | Size |
|--------|------|
| 556050 | 23G |
| 656050 | 25G |
| 856050 | 27G |



ICL Forceps

651072


| SKU | Size |
|--------|------|
| 651072 | 21G |



IOL Cutting Forceps

651075

| SKU | Size |
|--------|------|
| 651075 | 23G |



ICL Forceps Serrated Without Hole

651078

| SKU | Size |
|--------|------|
| 651078 | 21G |



MICRO INSTRUMENTS

Micro Forceps

ICL Forceps Serrated With Hole

651079

| SKU | Size |
|--------|------|
| 651079 | 21G |



Corneal Endothelium Forceps

651092

| SKU | Size |
|--------|------|
| 651092 | 23G |



MICRO INSTRUMENTS

Micro Scissors

Vitreoretinal Scissors, Straight Blunt Tips

556001

| SKU | Size |
|--------|------|
| 556001 | 23G |
| 656001 | 25G |
| 856001 | 27G |



Vitreoretinal Scissors, Straight Pointed Tips

556003

| SKU | Size |
|--------|------|
| 556003 | 23G |
| 656003 | 25G |
| 856003 | 27G |



Vitreoretinal Scissors, Wide Straight Pointed Tips

556004

| SKU | Size |
|--------|------|
| 556004 | 23G |
| 656004 | 25G |
| 856004 | 27G |

MOQ



Vitreoretinal Scissors, Curved

556006

| SKU | Size |
|--------|------|
| 556006 | 23G |
| 656006 | 25G |
| 856006 | 27G |



Vitreoretinal Scissors, Angled 45°

556008

| SKU | Size |
|--------|------|
| 556008 | 23G |
| 656008 | 25G |
| 856008 | 27G |



MICRO INSTRUMENTS

Micro Scissors

Vitreoretinal Scissors, Vertical 80°

556010

| SKU | Size |
|--------|------|
| 556010 | 23G |
| 556011 | 25G |



Vitreoretinal Scissors, Vertical 40°

556013

| SKU | Size |
|--------|------|
| 556013 | 23G |
| 556014 | 25G |



Vitreoretinal Scissors, Curved Shafts

556017

| SKU | Size |
|--------|------|
| 556017 | 23G |
| 656017 | 25G |
| 856017 | 27G |



Vitreoretinal Scissors, Long Curved Shafts

556018

| SKU | Size |
|--------|------|
| 556018 | 23G |
| 656018 | 25G |
| 856018 | 27G |



IOL Cutting Scissors 19G

653059

| SKU | Size |
|--------|------|
| 653059 | 19G |



MICRO INSTRUMENTS

Micro Capsulorhexis

Capsulorhexis Forceps, With Cystotome, Long Ending, 1.6 mm Platforms, Curved, 23G

556023

| SKU | Size |
|--------|------|
| 556023 | 23G |
| 656023 | 25G |
| 856023 | 27G |



Capsulorhexis Forceps, Without Cystotome, Curved Shaft

556040

| SKU | Size |
|--------|------|
| 556040 | 23G |
| 656040 | 25G |
| 856040 | 27G |



Capsulorhexis Forceps, With Cystotome, Short Ending, Curved Shaft

556041

| SKU | Size |
|--------|------|
| 556041 | 23G |
| 656041 | 25G |
| 856041 | 27G |



Trabeculectomy Punch

554023

| SKU | Length Tip |
|--------|------------|
| 554020 | 0.75 mm |
| 554021 | 1.0 mm |
| 554022 | 1.5 mm |
| 554023 | 2.0 mm |



Retinal Membrane Pick, 0.5mm Tip

554045

| SKU | Length Tip |
|--------|------------|
| 554045 | 0.5 mm |



Retinal Membrane Pick, Rounded Tip Spatula

554046

| SKU | Length Tip |
|--------|------------|
| 554046 | 0.5 mm |



Retinal Membrane Pick, Tip Angled 110°

554047

| SKU | Length Tip |
|--------|------------|
| 554047 | 0.5 mm |



Retinal Membrane Pick, 0.5mm Tips

554048

| SKU | Length Tip |
|--------|------------|
| 554048 | 0.5 mm |






SPECULUMS

SPECULUMS

Wire Speculum


Naqahara Wire Speculum,
13.5mm Open Blades

555000



Wire Speculum,
9mm Open Blades

555005



Kratz Barraquer Wire Speculum,
13.5mm Open Blades

555001




Wire Speculum,
14mm Open Blades

555006



Wire Speculum,
13.5mm Open Blades

555002




Wire Speculum,
14mm Open Blades

555007




Wire Speculum,
13.5mm Open Blades

555003




Wire Speculum,
9mm Open Blades

555008



Wire Speculum,
6mm & 9mm Open Blades

555004



Wire Speculum,
11mm Open Blades

555009



SPECULUMS

Wire Speculum

Jaffe Lid Retractors x2,
16mm Blades

555010



Barraquer Solid Speculum,
14.5mm Open Blades

555015



Jaffe Lid Retractors x2,
12mm Blades

555011



Wire Speculum,
10mm Open Blades

555016



Sprung Speculum Paediatric,
6mm Open Blades, L 26.5mm

555012



Wire Speculum,
14mm Open Blades

555017




Sprung Speculum Neonate,
4mm Open Blades, L 25mm

555013



Wire Speculum,
13mm Open Blades

555018




Barraquer Solid Speculum,
14mm Open Blades

555014



Wire Speculum,
14mm Open Blades

555019



SPECULUMS

Wire Speculum


Pierce Type Adjustable Speculum,
16mm Open Blades

555020



Paediatric Wire Speculum,
9mm Open Blades

655008




SPECULUMS

Adjustable Speculum

Adjustable Speculum, Swivel Open Blades

555022



Lieberman Adjustable Speculum, V-Shaped Wire

555023



Reversible Speculum, Lieberman Speculum, K-Wire, 13mm Open Blades, 45° Angled Shaft

555024



Kershner Eye Speculum, 13.5mm Open Blades

555025



Adjustable Speculum, Steinert Lasik Speculum, 14mm Inner Attached Round Solid Blades

555026



SPECULUMS

Adjustable Speculum

Adjustable Speculum, Speculum,
13mm Outer Attached Solid Square Blades

555027



Adjustable Speculum, Speculum,
14mm Round Solid Blades

555028



Adjustable Speculum, Lieberman Speculum,
12.5mm Inner Attached Solid Square Blades

555029




Adjustable Speculum, Ginsberg Eye Speculum,
13mm Solid Blades

555030



Lindstrom-Chu Aspirating Speculum

555031



SPECULUMS

Adjustable Speculum

Adjustable Speculum, Lieberman Speculum,
13mm Open Blades

555033



Adjustable Speculum, Lieberman Speculum,
14mm V Type Solid Blades, Wire Type

555034



Speculum Adjustable, V-Wire Angled

655023



Speculum Lieberman Fine, Open Blades

655035



Speculum Adjustable Full Squared Blades

655037



SPECULUMS

Adjustable Speculum

Speculum Adjustable Fine Full Blades

655038




Speculum Lieberman Open Rounded Blades

655040



Speculum Adjustable Full Rounded Blades

655041



Solid Bladed Speculum Screw Lock

655042





CURETTES & RETRACTORS

CURETTES & RETRACTORS

Curette, Double Ended

554030



Curette


554031

| SKU | Diameter |
|--------|----------|
| 554031 | 1.5mm |
| 554032 | 2mm |
| 554033 | 2.5mm |



Arruga Globe Retractor
10.5mm Wide Retractor Blade & 2.8mm Wide Retractor Blade

554080



Desmarres Hook Retractor

554090

| SKU | Size |
|--------|--------|
| 554090 | Small |
| 554091 | Medium |
| 554092 | Large |



Lachrymal Sac Retractor,
Curved 0.3×10mm Teeth

555032





DEPRESSORS AND HOOKS

DEPRESSORS & HOOKS

2.0mm Tear Drop Shape Tip,
Scleral Depressor

554024



6mm Wide Bar,
Scleral Depressor

554025



Scleral Depressor,
Schocket Scleral Depressor

554026



Jameson Strabismus Hook

554069

| SKU | Hook |
|--------|-------|
| 554068 | 7mm |
| 554069 | 8.5mm |



Chavasse Hook, Helveston Muscle Hook,
12mm Hook

554070



Retinal Detachment Hook,
1.5mm Oval Hold

554072



DEPRESSORS & HOOKS

Gillis Hook

554086



Indentor

654024




Muscle Hook

654070

| SKU | Diameter |
|--------|---------------|
| 654070 | Short Ending |
| 654071 | Long Ending |
| 654072 | Medium Ending |





MISCELLANEOUS



MISCELLANEOUS

Bulldog Clips, Flat Handle, 40mm

552038



Blade Breaker, Straight,
Round Knurled Handle With End Lock

552039



Blade Breaker, Curved,
Round Knurled Handle With End Lock

552040

MOQ



Blade Handle

552041



Freer Lacrimal Chisel,
Single Ends With 5.8×12mm

554049



MISCELLANEOUS

Fine Thornton Fixation Ring, Swivel Head

554074

| SKU | Diameter | Cut Out |
|--------|----------|---------|
| 554074 | 16mm | 14mm |
| 554075 | 13mm | 11mm |



Cautery, 5x7mm Oval Ball Sharp Tip With 6mm Long Point


554081



Foreign Body Remover

651077






BESTSELLERS

BESTSELLERS


Toothed Forceps, Straight Shafts With 5mm Tying Platforms,
Flat Handle, 85mm, 0.12mm, 1×2 Teeth

551001



Tying Forceps, Curved Shafts With 6mm Tying Platform,
Flat Handle, 85mm

551175




Colibri Toothed Forceps, Curved Shafts, Serrated Handle,
0.12mm, 1×2 Teeth

551187




Jeweler Style Forceps, Angled 45°, 5#, Extra Delicate,
0.05×0.01mm Tip

551222



Moorfields Forceps, Blunt Tipped Forceps

551235



Vannas Scissors, Curved, Sharp Pointed Tips,
Tip To Pivot Length 10mm

553014



Haefliger Phaco Chopper (Chopper & Manipulator in One)

554064



Vitreoretinal Forceps, Short End-Gripping Forceps (ILM), 23G

556026



BESTSELLERS

Mendez Degree Gauge,
Calibrated Every 10° From 0° to 180°,
External Diameter 19.0mm, IØ 12.0mm

558002



Capsulorhexis Forceps, Curved Thin Shafts,
Round Knurled Handle

651141



Iris Scissors Hasa Optix,
Straight, Sharp

653041



Lenticule Separator/Elevator

654108



LASIK Flap Lifter/Sinskey

654124



Speculum Lieberman Fine, Open Blades

655035






Speculum Adjustable Full Rounded Blades

655041



Contact us for e-catalogues and discover HASA OPTIX surgical instrument solutions




info@hasaoptix.com



hasaoptix.com



+3225246388



HASA OPTIX
Belgium


Join us to simplify
ophthalmic surgery.



Contact us to find your nearest distributor

 info@hasaoptix.com

 hasaoptix.com

 +32.2.524.63.88

523 Avenue Louise
1050, Brussels, Belgium

Follow us    

