



OPTIMISE COST AND QUALITY OF **CORNEA** SURGERY

Our innovative surgical instrument solution provides total peace of mind, confidence and satisfaction.



HASA OPTIX
Belgium

Build your own sets for **CORNEA** surgeries



PREMIUM QUALITY

Premium quality disposable recyclable instruments offering the constant quality with increased convenience and efficiency.

COST CONTROL

HASA OPTIX instruments provide the ideal procurement solution by aiding the precise calculation of surgery costs.

SIMPLIFY

Our pre-packed pouches help reduce the number of SKUs, contributing towards a simpler, more efficient repository solution.

PEACE OF MIND


Our innovative solution provides total peace of mind of always having accessible right instruments.

INFECTION CONTROL

Our single use and pre-sterilized instruments ensure all surgeons have peace of mind regards patient safety.

ENVIRONMENTALLY FRIENDLY

Supplied in recyclable pre-packed pouches relevant to each ophthalmic procedure, helping to reduce waste and environmental impact.




HASARECYCLE
EXCELLENCE PROGRAM




At HASA OPTIX, we have decided to create a fully integrated re-use circle so that the stainless steel, our instruments are made of, can have a second life as non-medical objects.


Anatomy Forceps



Anatomy Scissors



Needle Holders



OVERALL LENGTH

TIPS

- Fine / Micro**
for finest sutures
- Standard**
for general purpose
- Heavy**
for larger sutures

SHAFTS

- Straight***
- Curved***

*Choice depends on users preference although curved shafts may give better visualisation of the wound.

HANDLES


- Flat**
for enhanced manual grip and control of the tip (significant hand, wrist, or forearm movement may be required)
- Round**
enhances micro control of the tips with minimum movement by rolling the handle between the fingers
- Lock Option**
used to lock the jaws in the closed position, securing the needle during suturing

HASA OPTIX'S BESTSELLERS


Speculums



655035




555022




555025



655041



555023



655040

- 655035 Lieberman Speculum
- 555022 Adjustable Speculum (Swivel)
- 555025 Closed Wire Adjustable Blades
- 655041 Solid Round Adjustable Blades
- 555023 Lieberman V-wire Adult Adjustable
- 655040 Open Wire Round Adjustable Blades

Calipers




558009



554038

- 558009 Castroviejo Caliper
- 554036 Alcohol Well 7.0mm
- 554037 Alcohol Well 8.0mm
- 554038 Alcohol Well 9.0mm

Forceps



551119



551187



551175



551169



551122



551001



651082

- 551119 Polack Corneal Suturing Forceps
- 551187 Colibri Toothed Forceps 0.12mm Serrated Handle
- 551175 Tying Forceps Curved Jaws Flat Handle
- 551169 Tying Forceps Platform
- 551122 Notched Forceps Straight
- 551001 Toothed Forceps 0.12mm Flat Handle
- 651082 Corneal Forceps 0.20mm

Manipulators



654123

- 654123 Doubled-ended Reversed Sinskey/Sinskey

Micro Instruments




651075



651079



556028



653059

- 651075 IOL Cutting Forceps
- 651079 ICL Forceps Serrated with Hole
- 556028 VR Round Short Flat Forceps 23G
- 653059 IOL Cutter 19G

Scissors



553009



553010




553014



553036

- 553009 Castroviejo Corneal Scissors 9mm (L)
- 553010 Castroviejo Corneal Scissors 9mm (R)
- 553014 Vannas Scissors Curved
- 553036 Vannas Scissors Straight

Needle Holders




552027



552022



552025



552005

- 552027 6mm Curved Jaws Round Knurled Handle
- 552022 10mm Curved Jaws Round Knurled Handle with Lock
- 552025 9mm Curved Jaws Round Knurled Handle
- 552005 6mm Straight Jaws Round Knurled Handle

Cataract

Cornea

Refractive

Glaucoma

Oculoplastic

Vitreoretinal

SIMPLIFYING OPHTHALMIC SURGERY



Contact us to find your nearest distributor


@ info@hasaoptix.com

www.hasaoptix.com

+32.2.524.63.88

523 Avenue Louise
1050, Brussels, Belgium

Follow us



HASARECYCLE
EXCELLENCE PROGRAM

