A genuine desire to empower all those involved in ophthalmic surgery


HASA OPTIX holds a genuine desire to empower all those involved in ophthalmic surgery with precision performance, procurement reliability, and patient quality of care always at the forefront of our approach.

We design and manufacture premium, single-use, and fully recyclable ophthalmic instruments for both anterior and posterior segment surgeries.

Made of stainless steel, our exceptional precision instruments are supplied in recyclable pre-packed pouches relevant to each ophthalmic procedure. This not only ensures all surgeons have peace of mind regards sterilisation and reliability, but also provides a cost controlled, simplified procurement and repository solution.

We are passionate about eliminating cross contamination and reducing risk management whilst increasing efficiency, convenience, and reliability.

Our innovative solution provides total peace of mind, confidence, and satisfaction.





Certified B Corporation

MOST SUSTAINABLE COMPANY 2025

What Sets Us Apart?

Quality improvement and cost-savings are of great importance in healthcare in general and in ophthalmic surgery. Hospitals, more than ever, must face the challenges of working profitably and obtaining optimum quality at the same time. In a changing world, we believe there are new ways that will allow to facilitate work of all stakeholders involved in ophthalmic surgery, allowing a higher pharmaco economic value.

Premium quality

Premium quality single use recyclable instruments offering the same constant quality as reusable instruments but with increased convenience and efficiency.

Simplicity

Our pre-packed pouches help reduce the number of SKUs, contributing towards a simpler, more efficient repository solution.

Cost control

HASA OPTIX instruments provide the ideal procurement solution by aiding the precise calculation of surgery costs.

Peace of mind

Our innovative solution provides total peace of mind ensuring the correct surgical instrumentals are easily accessible.

Infection control

Our single use and pre-sterilised instruments ensure all surgeons can be confident with their patient safety of care.

Environmentally friendly


Supplied in recyclable pre-packed pouches relevant to each ophthalmic procedure, helping to reduce waste and environmental impact.

Important Information:
Individual instruments are supplied in a box of 10 pcs, unless otherwise specified.
Instrument sets are supplied in a box of 10 sets, unless otherwise specified.
When MOQ is indicated, it signifies that the Minimum Order Quantity for that SKU is 1000 pieces.

Contents

Forceps 4
Needle Holders 40
Manipulators 46
Choppers 56
Spatulas 62
Double-Ended Instruments 64
Probes 72
Dilators & Loops 74
Markers 76
Gauges 78
Calipers 80
Scissors 82
Vitreoretinal Forceps 96
Vitreoretinal Scissors 102
Micro Instruments 106
Speculums 110
Curettes 118
Retractors 120
Depressors 122
Hooks 124
Miscellaneous 126
Bestsellers 130

Forceps



Capsulorhexis Forceps

651000

Capsulorhexis Forceps
Curved Thin Shaft, Flat Cross-Action Handle



651001

Capsulorhexis Forceps
45° Angled Thin Shaft, Flat Cross-Action Handle



651002

Capsulorhexis Forceps
Curved Thin Shaft,
Cross-Action, Round Knurled Handle



Capsulorhexis Forceps

651003

Capsulorhexis Forceps
45° Angled Thin, Cross-Action Shaft,
Round Knurled Handle



651141

Capsulorhexis Forceps
Curved Thin Shaft, Round Knurled Handle



651143


Capsulorhexis Forceps
Curved Thin Shaft,
Cross-Action Ripped Handle



Capsulorhexis Forceps

651147

Capsulorhexis Forceps
45° Angled Thin Shaft, Round Knurled Handle



651149

Capsulorhexis Forceps
45° Angled Thin Shaft,
Cross-Action Ripped Handle



651150

Capsulorhexis Forceps
45° Angled Thin Shaft, Flat Handle

SKU	Length
651150	105 mm
651151	85 mm
651152	115 mm

MOQ



651153

Capsulorhexis Forceps
Curved Thin Shaft, Flat Handle

SKU	Length
651153	85 mm
651154	105 mm
651155	115 mm

MOQ



Fixation Forceps

551000

Toothed Forceps
Straight Shafts, 5 mm Tying Platforms,
Flat Handle, 85 mm Length, 1×2 Teeth

SKU	Size
551000	0.10 mm
551001	0.12 mm
551002	0.15 mm
551003	0.30 mm
551004	0.40 mm
551005	0.50 mm



551006

Toothed Forceps
Straight Shaft, 5 mm Tying Platforms,
Flat Handle, 105 mm Length, 1×2 Teeth

SKU	Size
551006	0.10 mm
551007	0.12 mm
551008	0.15 mm
551009	0.30 mm
551010	0.40 mm
551011	0.50 mm



551012

Toothed Forceps
Straight Shaft, 5 mm Tying Platforms,
Flat Handle, 115 mm Length, 1×2 Teeth

SKU	Size
551012	0.10 mm
551013	0.12 mm
551014	0.15 mm
551015	0.30 mm
551016	0.40 mm
551017	0.50 mm



Forceps

Fixation Forceps

551018

Toothed Forceps
Straight Shaft, 5 mm Tying Platforms,
Round Knurled Handle, 1×2 Teeth

SKU	Size
551018	0.10 mm
551019	0.12 mm
551020	0.15 mm
551021	0.30 mm
551022	0.40 mm
551023	0.50 mm



551024

Toothed Forceps
Straight Shaft, 5 mm Tying Platforms,
Castroviejo Handle, 1×2 Teeth


SKU	Size
551024	0.10 mm
551025	0.12 mm
551026	0.15 mm
551027	0.30 mm
551028	0.40 mm
551029	0.50 mm



551037

Curved Toothed Forceps
Curved Shaft, 5 mm Tying Platforms,
Flat Handle, 1×2 Teeth

SKU	Size
551037	0.10 mm
551038	0.12 mm
551039	0.15 mm
551040	0.30 mm
551041	0.40 mm
551042	0.50 mm



Fixation Forceps

551043

Colibri Conjunctiva Toothed Forceps
Curved Shaft, 5 mm Tying Platforms,
Flat Handle, 1×2 Teeth

SKU	Size
551043	0.10 mm
551044	0.12 mm
551045	0.15 mm
551046	0.30 mm
551047	0.40 mm
551048	0.50 mm



551049

Colibri Toothed Forceps
Curved Shaft, 5 mm Tying Platforms,
Flat Handle, 1×2 Teeth

SKU	Size
551049	0.10 mm
551050	0.12 mm
551051	0.15 mm
551052	0.30 mm
551053	0.40 mm
551054	0.50 mm



551055

Colibri Toothed Forceps
Curved Shaft, 5 mm Tying Platforms,
Flat Handle, 1×2 Teeth

SKU	Size
551055	0.10 mm
551056	0.12 mm
551057	0.15 mm
551058	0.30 mm
551059	0.40 mm
551060	0.50 mm



Fixation Forceps

551061

Curved Toothed Forceps
Curved Shaft, 5 mm Tying Platforms,
Round Knurled Handle, 1×2 Teeth

SKU	Size
551061	0.10 mm
551062	0.12 mm
551063	0.15 mm
551064	0.30 mm
551065	0.40 mm
551066	0.50 mm



551067

Curved Toothed Forceps
Curved Shaft, 5 mm Tying Platforms,
Castroviejo Handle, 1×2 Teeth

SKU	Size
551067	0.10 mm
551068	0.12 mm
551069	0.15 mm
551070	0.30 mm
551071	0.40 mm
551072	0.50 mm



551080

Kelman-McPherson Forceps
Angled Shaft, 5 mm Tying Platforms,
Flat Handle, 85 mm Length, 1×2 Teeth

SKU	Size
551080	0.10 mm
551081	0.12 mm
551082	0.15 mm
551083	0.30 mm
551084	0.40 mm
551085	0.50 mm



Fixation Forceps

551086

Kelman-McPherson Forceps
Angled Shaft, 5 mm Tying Platforms,
Flat Handle, 105 mm Length, 1×2 Teeth

SKU	Size
551086	0.10 mm
551087	0.12 mm
551088	0.15 mm
551089	0.30 mm
551090	0.40 mm
551091	0.50 mm



551092

Kelman-McPherson Forceps
Angled Shaft, 5 mm Tying Platforms,
Flat Handle, 115 mm Length, 1×2 Teeth

SKU	Size
551092	0.10 mm
551093	0.12 mm
551094	0.15 mm
551095	0.30 mm
551096	0.40 mm
551097	0.50 mm



551098

Kelman-McPherson Forceps
Angled Shaft, 5 mm Tying Platforms,
Round Knurled Handle, 1×2 Teeth

SKU	Size
551098	0.10 mm
551099	0.12 mm
551100	0.15 mm
551101	0.30 mm
551102	0.40 mm
551103	0.50 mm



Forceps

Fixation Forceps

551104

Kelman-McPherson Forceps
Angled Shaft, 5 mm Tying Platforms,
Castroviejo Handle, 1×2 Teeth

SKU	Size
551104	0.10 mm
551105	0.12 mm
551106	0.15 mm
551107	0.30 mm
551108	0.40 mm
551109	0.50 mm



551119

Polack Corneal Suturing Forceps
Curved Shaft, Flat Handle,
0.12 mm Size, 1×2 Teeth



551122

Notched Forceps
Straight Shaft, Tying Platforms,
Flat Handle, 0.25 mm Size


SKU	Length
551122	85 mm
551123	105 mm
551124	115 mm



Fixation Forceps

551125

Notched Forceps
Straight Shaft, Tying Platforms,
Round Knurled Handle, 0.25 mm Size, Notched



551126

Notched Forceps
Straight Shaft, Tying Platforms,
Castroviejo Handle, 0.25 mm Size, Notched



551128

Notched Forceps
Curved Shaft, Tying Platforms,
Flat Handle, 0.25 mm Size, Notched

SKU	Length
551128	85 mm
551129	105 mm
551130	115 mm



Forceps

Fixation Forceps

551131

Pierse-Colibri Corneal Forceps
Curved Shaft, Tying Platforms,
Flat Serrated Handle, 0.25 mm Size, Notched



551132

Notched Forceps
Curved Shaft, Tying Platforms, Round Handle,
0.25 mm Size, Notched



551135

Notched Forceps
45° Angled Straight Shaft, Tying Platforms,
Flat Handle, 0.25 mm Size, Notched

SKU	Length
551135	85 mm
551136	105 mm
551137	115 mm



Fixation Forceps


551138

Notched Forceps
45° Angled Straight Shaft, Tying Platforms,
Round Handle, 0.25 mm Size, Notched



551139

Notched Forceps
45° Angled Straight Shaft, Tying Platforms,
Castroviejo Handle, 0.25 mm Size, Notched



551141

Hofmann-Polack Corneal Suturing Forceps
45° Angled Straight Shaft, Flat Handle,
0.25 mm Size, Notched



Fixation Forceps

551187

Colibri Toothed Forceps
Curved Shaft, 0.12 mm, 1x2 Teeth



651080

Paufique Forceps
Tying Platforms, Oblique Teeth 0.7 mm



651081

Corneal Forceps
Straight Shaft, Tying Platforms,
Dipped Serrated Handle

SKU	Teeth
651081	Oblique 0.12 mm
651082	Oblique 0.20 mm



Fixation Forceps

651085

Corneal Forceps
Curved Shaft, Tying Platforms,
Dipped Serrated Handle, Oblique Teeth 0.12 mm



Forceps

Tissue Forceps


551120

Adson Forceps
Straight Shaft, No Platform,
Swallow-Tail Handle, 0.40 mm Size, 1×2 Teeth



551235

Moorfields Forceps
Blunt Tipped, For Manipulating the Conjunctiva



551240

Gillis Tissue Forceps

SKU	Size
551240	0.50 mm
551241	0.40 mm
551242	0.30 mm



Tissue Forceps

651120

Adson Forceps Non-Toothed
Serrated Jaws



Forceps

Tying Forceps

551169

Tying Forceps
Straight Shaft, 4.5 mm Tying Platform,
Flat Handle

SKU	Size
551169	85 mm
551170	105 mm
551171	115 mm



551172

Tying Forceps
Straight Shaft, 4.5 mm Tying Platform,
Round Knurled Handle



551175

Tying Forceps
Curved Shaft, 6 mm Tying Platform,
Flat Handle

SKU	Size
551175	85 mm
551176	105 mm
551177	115 mm



Tying Forceps

551178

Tying Forceps
Curved Shaft, 6 mm Tying Platform,
Round Knurled Handle



551179

Tying Forceps
Curved Shaft, 6 mm Tying Platform,
Castroviejo Handle



551181

Tying Forceps
45° Angled Shaft, 4.5 mm Tying Platform,
Flat Handle


SKU	Length
551181	85 mm
551182	105 mm
551183	115 mm



Tying Forceps


551184

Tying Forceps
45° Angled Shaft, 4.5 mm Tying Platform,
Round Knurled Handle



551185

Tying Forceps
45° Angled Shaft, 4.5 mm Tying Platform,
Castroviejo Handle



651187

McPherson Tying Forceps
7.5 mm Tying Platforms



Lens Forceps


551188

Lens Forceps
Curved Shaft, Cross-Action Handle,
90° Angled, 1 mm Tips



551189

Lens Forceps
Curved Shaft, Cross-Action Handle,
1 mm Ring Tips



551190

Lens Forceps
Angled Shaft, Flat Handle, 1 mm Round Tips



Forceps

Lens Forceps

551191

Lens Forceps
45° Angled Shaft, Cross-Action Handle,
Smooth Flat Jaws



551192


Kelman-McPherson Lens Forceps
Angled Shaft, Castroviejo Handle,
Smooth Jaws, 105 mm



551194

Kelman-McPherson Lens Forceps
Cross-Action Handle, Smooth, Curved Jaws

SKU	Jaws
551194	With Lock
551195	Without Lock



Lens Forceps

551199

Kelman-McPherson Lens Forceps
Angled Shaft, Flat Handle



651071

ICL Cartridge Loading Forceps
Cross-Action Shaft



651073

ICL Cartridge Loading Forceps
Flat Handle



Forceps

Lens Forceps

651076

IOL Forceps



651086


Dennis Forceps
Straight Shaft, Flat Handle

MOQ



651160

Lenticule Extraction Forceps



Lens Forceps

651161

Lenticule Removal Forceps
Angled Shaft, Castroviejo Handle,
Serrated Jaws



LASIK Forceps

551211

Perone LASIK Flap Forceps
2 mm Rings Tips

551212

Perone LASIK Flap Forceps
2 mm Long Tips



551225

LASIK Retreatment Forceps
1×2 Teeth, 0.50 mm Size

551230


LASIK Multipurpose Forceps
Curved Jaws



Pre-Chopper Forceps

551200

Akahoshi Prechopper Forceps
Flat Handle, Straight, Cross-Action Shaft,
Maximum Tip Opening 3 mm



551201

Akahoshi Prechopper Forceps
Round Handle, Straight, Cross-Action Shaft



551205

Akahoshi Prechopper Forceps
Cross-Action Handle, 11 mm Angled Long Shaft



Pre-Chopper Forceps


551206

Akahoshi Prechopper Forceps
Cross-Action Handle, 10 mm Angled Shaft



551207

Akahoshi Prechopper Forceps
Cross-Action Handle, 12 mm Angled Shaft



551208

Akahoshi Prechopper Forceps
Cross-Action Handle, 12 mm Curved Shaft



Chalazion Forceps

551213

Chalazion Forceps
Serrated, Locking Thumb Screw Handle,
Round 8 mm Plate



551214

Chalazion Forceps
Serrated, Locking Thumb Screw Handle,
Oval 8×10 mm Plate



551215

Chalazion Forceps
Serrated, Locking Thumb Screw Handle,
Round 15 mm Plate



551216

Chalazion Forceps
Serrated, Locking Thumb Screw Handle,
Round 12 mm Plate



Chalazion Forceps

551217

Chalazion Forceps
Serrated, Locking Thumb Screw Handle,
Oval 12×14 mm Plate



551218

Ciliary Forceps
Flat Handle



551236

Chalazion Forceps
Round 11 mm Plate



551237

Chalazion Forceps
Oval 21×12 mm Plate



General Forceps


551219

Scleral Plug Forceps
Cross-Action Shaft, Swallow-Tail Handle



551220

Jeweler Style Forceps
Straight Shaft, 0.05x0.01 mm Tip



551222

Jeweler Style Forceps
45° Angled Shaft, 0.05x0.01 mm Tip,
Extra Delicate



General Forceps

551223

Jeweler Style Forceps
Curved Shaft, 0.07×0.03 mm Tip



551224

Landers Vitrectomy Lens Forceps
45° Angled Shaft, Curved Grooved Blades



552034

Mosquito Forceps
Straight Shaft, 24 mm Serrated Jaws



General Forceps

552035

Mosquito Forceps
Curved Shaft, 24 mm Serrated Jaws



552042

Towel Clamps
Flat Shaft



651070

Meibomian Forceps



General Forceps

651090

Florakis Endothial Forceps

MOQ



651091

Iris Claw Forceps



651156

Mermoud Forceps



General Forceps

651162

Plain Forceps



Micro Capsulorhexis Forceps

556023

Capsulorhexis Forceps
With Cystotome, Curved Shaft, Long Ending

SKU	Gauge
556023	23G
656023	25G
856023	27G





556040

Capsulorhexis Forceps
Without Cystotome, Curved Shaft

SKU	Gauge
556040	23G
656040	25G
856040	27G





556041


Capsulorhexis Forceps
With Cystotome, Curved Shaft, Short Ending

SKU	Gauge
556041	23G
656041	25G
856041	27G





Needle Holders



Needle Holders

Needle Holders

552001

Needle Holder
Straight Jaws, 9 mm Jaw Size, Flat Handle

SKU	Lock
552000	Yes
552001	No



552003

Needle Holder
Straight Jaws, 11 mm Jaw Size, Flat Handle

SKU	Lock
552002	Yes
552003	No



552005

Needle Holder
Straight Jaws, 6 mm Jaw Size,
Round Knurled Handle

SKU	Lock
552004	Yes
552005	No



552007

Needle Holder
Straight Jaws, 9 mm Jaw Size, Flat Handle

SKU	Lock
552006	Yes
552007	No



Needle Holders

552009

Needle Holder
Curved Jaws, 11 mm Jaw Size,
Round Knurled Handle

SKU	Lock
552008	Yes
552009	No



552011

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Round Knurled Handle

SKU	Lock
552010	Yes
552011	No



552013

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Flat Handle

SKU	Lock
552012	Yes
552013	No



552015

Needle Holder
Curved Jaws, 11 mm Jaw Size,
Flat Handle

SKU	Lock
552014	Yes
552015	No



Needle Holders

552017

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Flat Handle

SKU	Lock
552016	Yes
552017	No



552018

Needle Holder
Curved Jaws, 9.0 mm Jaw Size,
Round Knurled Handle

SKU	Lock
552018	Yes
552019	No



552021

Needle Holder
Curved Jaws, 9 mm Jaw,
Round Knurled Handle

SKU	Lock
552020	Yes
552021	No



552023

Needle Holder
Curved Jaws, 10 mm Jaw Size,
Round Knurled Handle

SKU	Lock
552022	Yes
552023	No



Needle Holders

552025

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Round Knurled Handle


SKU	Lock
552024	Yes
552025	No



552027

Needle Holder
Curved Jaws, 6 mm Jaw Size,
Round Knurled Handle

SKU	Lock
552026	Yes
552027	No



552029

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Flat Handle

SKU	Lock
552028	Yes
552029	No



Needle Holders

Needle Holders

552030

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Round Knurled Handle, No Lock



552032

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Flat Handle, No Lock



652021

Needle Holder
Curved Jaws, 9 mm Jaw Size,
Round Knurled Handle, No Lock



Manipulators




Manipulators

Single-Ended Manipulators

554000

Kuglen Iris Hook
K-Shaped Tip


SKU	Tip
554000	Angled
554001	Straight



554002


Jaffe-Knolle Iris Hook
Iris Hook

SKU	Tip
554002	Angled
554003	Straight



554004


Pupil Snapper Hook
Angled Tip, Smoothly Polished Notch, Sharp Inner Edge



554005

Lens Manipulator
T-Shaped Tip

SKU	Tip
554005	Straight
554006	Angled




Single-Ended Manipulators

554007

Lens Manipulator
Y Nucleus Rotator, Angled


SKU	Shape
554007	Vertical
554087	Horizontal



554008


Lens Manipulator
V-Shaped Tip

SKU	Tip
554008	Angled
554009	Straight



554010

Lens Manipulator
Angled Tip, 0.7 mm Tip Size



Single-Ended Manipulators


554011

Lens Manipulator
Curved Tip, 1.1 mm Tip Size




554012

Nucleus Manipulator
Drysdale



554013


Nucleus Mushroom Manipulator



Single-Ended Manipulators


554014

Nucleus Manipulator
Connor Wand, Straight Tip




554015

Nucleus Manipulator
Connor Wand, Angled Tip, Round



554016

Nucleus Manipulator
Connor Wand, Angled Tip, Oval



Single-Ended Manipulators

554018


Iris Hook
H Tip, Angled Tip



554027


Repositor
Semi-Sharp Tip

SKU	Tip Size
554027	0.7 mm
554028	0.5 mm



554053


Kansas Nucleus Bisector



Single-Ended Manipulators

554071

Intra-LASIK Flap Separator
Round, Sharp Tipped Separator



554504


LASIK Flap Lifter/Elevator Binder



654073

ICL Manipulator

SKU	Shape	Tip Diameter
654073	Angled	1 mm
654074	Curved	1 mm
654075	Curved	0.6 mm
654076	Angled	0.6 mm




Manipulators

Single-Ended Manipulators


654088

Lund Iris Hook
45° Angled




654097

Manipulator
Lester Angled



654098


Sinskey Hook
Straight



Single-Ended Manipulators


654107

Lebuisson Manipulator




654112

Manipulator Cylinder



654113


Manipulator Govan



Single-Ended Manipulators

654119


Manipulator Clayman Spatula
Flat



654120


Manipulator Sinskey
0.2 mm Tip Size

SKU	Details
654120	0.7 mm Altitude
654121	1.1 mm Altitude



654128

Multifunctional Manipulator



Choppers



Choppers

Single-Ended Choppers

554054


Phaco Chopper

SKU	Hand
554054	Left Hand
554055	Right Hand




554056

Phaco Chopper
“M” Hook



554057

Phaco Chopper
Ambidextrous Nucleus Chopper



Single-Ended Choppers


554058

Phaco Chopper
Angled Nucleus Chopper




554059

Phaco Chopper
Curved Chopper,
2 Faceted Chopping Surfaces



554060

Phaco Chopper
Straight Blade, 1.8 mm Tip



Single-Ended Choppers

554061

Phaco Chopper
Curved Blade, 1.8 mm Tip



554062

Phaco Chopper
Angled Shaft, 1.6 mm Tip

SKU	Hand
554062	Left Hand
554063	Right Hand



554064


Haefliger Phaco Chopper
Chopper & Manipulator in One



Single-Ended Choppers


554065

Phaco Chopper
Special Inner Blade Offset 45°




554066

Kuglen Iris Hook
“S”, 0.6 mm Hook Tip



554067

Kuglen Iris Hook
Curved Hook Tip



Single-Ended Choppers

554088

Rosen Phaco Chopper



654114

Denis Nucleus Chopper



554085

Koch Stop & Chop Manipulator



Spatulas




Spatulas

Single-Ended Spatulas

554082

Corneal Epithelium Spatula




554083

Epi-Peeler



554084

Hoe Type Epi-Peeler
Square Shaped



554095

Koch Phaco Spatula



654129

Paufique Spatula




Double-Ended Instruments



Double-Ended Instruments


654086

Double-Ended Manipulator
Drysdale/Y Pusher




654087

Double-Ended Manipulator
Cylinder/Y Pusher



654109


Double-Ended Manipulator
Sinskey/Lester



Double-Ended Instruments


654118

Double-Ended Manipulator
Sweep/Y




654123

Double-Ended Manipulator
Reversed Sinskey/Sinskey



554501

Double-Ended Manipulator
Drysdale/Sinskey



Double-Ended Instruments

Double-Ended Instruments


554502

Double-Ended Manipulator
Fukasaku/Sinskey




554503

Double-Ended Manipulator
LASIK Flap Separator/Elevator



654085


Double-Ended Manipulator
Drysdale/H-Shaped Tip



Double-Ended Instruments


654089

Double-Ended Manipulator
Drysdale/Fukasaku Small Pupil Snapper




654090

Double-Ended Manipulator
Agarwal Phaco Chopper/
Fukasaku Small Pupil Snapper



654096

Double-Ended Manipulator
Drysdale/Koch Chopper



Double-Ended Instruments


654100

Double-Ended Manipulator
Drysdale/Curved Chopper




654104

Double-Ended Manipulator
Fukasaku Small Pupil Snapper/Cylinder



654106


Double-Ended Manipulator
Koch Chop/Chopper



Double-Ended Instruments


654108

Double-Ended Manipulator
Lenticule Separator/Elevator




654110

Double-Ended Manipulator
Lund Iris Hook/Phaco Chopper



654111

Double-Ended Manipulator
Mushroom/Chopper




Double-Ended Instruments

Double-Ended Instruments


654124

Double-Ended Manipulator
LASIK Flap Lifter/Sinskey



554029

Double-Ended Manipulator
Repositor



Probes



Probes

Probes

731

Lacrimal Probe, Set of 4 pieces
0.5/0.6 mm, 0.7/0.8 mm, 0.9/1 mm
& 1.1/1.2 mm Size



732

Lacrimal Probe
Pigtail Curved Probe



SKU	Diameter
732	6 mm
733	7 mm
734	8 mm
735	9 mm

Dilators & Loops



Dilators & Loops

Dilators & Loops

736

Punctal Dilator
Angled Nucleus Chopper

SKU	Taper	Tip Diameter
736	Sharp 16 mm	0.3 mm
737	Sharp 15 mm	0.4 mm
738	Sharp 15 mm	0.3 mm



739

Punctal Dilator
Sharp 17 mm Taper & Sharp 13 mm Taper, Double Ended



741

Double Ended Dilator
15 mm & 19 mm Taper



554034

Lens Loop

SKU	Loop
554034	5×7 mm
554035	3×8 mm



Markers



Markers

Markers

554036

LASEK Alcohol Well



SKU	Diameter
554036	7 mm
554037	8 mm
554038	9 mm
654039	9.5 mm

558000

Optical Zone, Target Centre
4 mm Diameter



MOQ

558001

CK LASIK Marker, Target Centre
9.5 mm Diameter



MOQ

558003

Double-Ended Scleral Marker
Aphakic 3.5/4 mm Diameter




558004

LASIK Corneal Marker
Two Blunt Fins, Separated at 45°Angle



Gauges



Gauges

558002

Mendez Degree Gauge
12 mm Diameter,
Calibrated Every 10 Degrees From 0° To 180°



Calipers



Calipers

558005

Braunstein Fixed Caliper
Aphakic 3.5/4 mm Diameter, Fine Tips Mark
Sclerotomy Sites for Vitreous Surgery



558006

Braunstein Fixed Caliper
Aphakic 3.5/4 mm Diameter, Fine Tips Mark
Sclerotomy Sites for Vitreous Surgery




558009

Castroviejo Caliper
Measure 0 mm To 20 mm Diameter,
Scale Reads On One Side



Scissors



Scissors

553000

Castroviejo Scissors
Curved, Blunt Tips,
14 mm Tip To Pivot Length



553001

Castroviejo Scissors
Straight, Blunt Tips,
14 mm Tip To Pivot Length



553002

Conjunctiva Scissors
Curved, Blunt Tips

SKU	Tip To Pivot Length
553002	11 mm
553003	12 mm
553004	16.5 mm



553005

Castroviejo Scissors
Curved, Blunt Tips

SKU	Blades	Tip To Pivot Length
553005	Large	16 mm
553006	Small	16.5 mm



Scissors

553007

Castroviejo Corneal Scissors
Blunt Tips, 13 mm Tip To Pivot Point



SKU	Blades
553007	Left, Lower
553008	Right, Lower

553009


Castroviejo Corneal Scissors
Blunt Tips, 11 mm Tip To Pivot Point



SKU	Blades
553009	Left, Lower
553010	Right, Lower

553011

Castroviejo Corneal Scissors
Curved, Blunt Tips



SKU	Tip To Pivot Length
553011	12 mm
553012	16 mm

Scissors

553013

Castroviejo Corneal Scissors
Angled, Lower Blades, Blunt Tips,
13 mm Tip To Pivot Point



553014

Vannas Scissors
Curved, Sharp Pointed Tips



SKU	Tip To Pivot Length
553014	10 mm
553015	16 mm

553016

Vannas Scissors
Straight, Sharp Pointed Tips,
10 mm Tip To Pivot Length



Scissors

553017

Vannas Scissors
Angled, Sharp Pointed Tips,
9 mm Tip To Pivot Point



553018

Iris Scissors
Curved, Sharp Tips,
6 mm Tip To Pivot Point, 11 mm Blade Length



553019

Iris Scissors
Curved, Sharp Tips,
13 mm Tip To Pivot Point, 13 mm Blade Length



Scissors

553020

Iris Scissors
Straight, Sharp Tips,
11 mm Blade Length, 6 mm Cutting Length



553021

Iris Scissors
Straight, Sharp Tips,
12 mm Blade Length



553022

Westcott Stitch Scissors
Curved, Sharp Pointed Tips

SKU	Tip To Pivot Length
553022	14 mm
553023	16 mm
553024	11 mm



Scissors

553025

Vannas Capsulotomy Scissors
Curved, Sharp Pointed Tips



SKU	Tip To Pivot Length
553025	11.5 mm
553026	15 mm

553027

Vannas Capsulotomy Scissors
Straight, Sharp Pointed Tips



SKU	Tip To Pivot Length
553027	13 mm
553028	16 mm

553029

Vannas Capsulotomy Scissors
Angled, Sharp Pointed Tips



SKU	Tip To Pivot Length
553029	12 mm
553030	14 mm

Scissors

553031

Tenotomy Scissors
Curved, Sharp Pointed Tips



SKU	Tip To Pivot Length	Overall Length
553031	16 mm	115 mm
553032	20 mm	125 mm

553034

Westcott Stitch Scissors
Curved, Sharp Pointed Tips,
16 mm Tip To Pivot Length



553035

Vannas Capsulotomy Scissors
Angled, Sharp Pointed Tips,
16 mm Tip To Pivot Length



Scissors

553036

Vannas Scissors
Straight, Sharp Pointed Tips,
10 mm Tip To Pivot Length



553038

Eye Scissors
Straight, Sharp Pointed Tips



553039

Eye Scissors
Straight, Blunt Ended Tips



SKU	Overall Length
553039	90 mm
553040	130 mm

Scissors

653034

Westcott Stitch Scissors
Curved, Sharp Tips, 21 mm Tip To Pivot Length
* for other variants of this SKU, please refer to pages 92-93



653041

Iris Scissors
Straight, Sharp Tips



653042

Iris Scissors
Curved, Sharp Tips



Scissors

653043

Iris Scissors
Curved, Blunt Tips



653044

Iris Scissors
Straight, Blunt Tips



653045

Westcott Stitch Scissors
Straight, Blunt Tips



Scissors

653055

Westcott Stitch Scissors
Curved, Blunt Tips,
22 mm Tip To Pivot Length



653056

Westcott Stitch Scissors
Straight, Sharp Tips
24 mm Tip To Pivot Length



653058

IOL Cutter



Scissors

653060

Noyes Iris Scissors
Straight, Blunt Tips



653061

Noyes Iris Scissors
Straight, Sharp Tips



653062

Noyes Iris Scissors
Curved, Blunt Tips



Scissors

653063

Noyes Iris Scissors
Curved, Sharp Tips



753013

Stevens Tenotomy Scissors
Straight, Short Blades, Blunt Tips
34 mm Tip to Pivot Length

MOQ



753035

Enucleation Scissors
Strong Curve, 40 mm Tip to Pivot Length

MOQ



Vitreoretinal Forceps



Vitreoretinal Forceps

Vitreoretinal Forceps

556026

Vitreoretinal Forceps
Short End-Gripping Forceps (ILM)

SKU	Gauge
556026	23G
656026	25G
856026	27G



556028

Vitreoretinal Forceps
Gripping Forceps, Smooth Rounded Tip

SKU	Gauge
556028	23G
656028	25G
856028	27G



556030

Vitreoretinal Forceps
Straight Shaft, Serrated Tips

SKU	Gauge
556030	23G
656030	25G
856030	27G




Vitreoretinal Forceps

556032

Vitreoretinal Forceps
Straight Shaft, 1×2 Teeth

SKU	Gauge
556032	23G
656032	25G
856032	27G



556035

Vitreoretinal Forceps
Membrane Peeler, Smooth, 30° Tips


SKU	Gauge
556035	23G
656035	25G
856035	27G



556037

Vitreoretinal Forceps
Membrane Peeler, Grooved, 30° Tips

SKU	Gauge
556037	23G
656037	25G
856037	27G



Vitreoretinal Forceps

556038

Vitreoretinal Forceps
Long End-Gripping Forceps

SKU	Gauge
556038	23G
656038	25G
856038	27G



556039

Vitreoretinal Forceps
End-Gripping Forceps, Round Tips

SKU	Gauge
556039	23G
656039	25G
856039	27G



556042

Vitreoretinal Forceps
Asymmetric End-Gripping

SKU	Gauge
556042	25G
856042	27G



Vitreoretinal Forceps

556043

Vitreoretinal Forceps
Pointed Short Forceps, Flat Tips


SKU	Gauge
556043	23G
656043	25G
856043	27G



556044

Vitreoretinal Forceps
Smooth Serrated Tips

SKU	Gauge
556044	23G
656044	25G
856044	27G




Vitreoretinal Forceps

556050

Vitreoretinal Forceps
Wide End-Gripping Forceps (ILM)

SKU	Gauge
556050	23G
656050	25G
856050	27G



651092

Corneal Endothelium Forceps
23G Gauge

MOQ



Vitreoretinal Scissors



Vitreoretinal Scissors

Vitreoretinal Scissors

556001

Vitreoretinal Scissors
Straight Blunt Tips

SKU	Gauge
556001	23G
656001	25G
856001	27G



556003

Vitreoretinal Scissors
Straight Pointed Tips

SKU	Gauge
556003	23G
656003	25G
856003	27G



556004

Vitreoretinal Scissors
Wide Straight Pointed Tips

SKU	Gauge	
556004	23G	
656004	25G	
856004	27G	MOQ



Vitreoretinal Scissors

556006

Vitreoretinal Scissors
Curved Tips

SKU	Gauge
556006	23G
656006	25G
856006	27G



556008

Vitreoretinal Scissors
Angled 45° Tips

SKU	Gauge
556008	23G
656008	25G
856008	27G



556010

Vitreoretinal Scissors
Vertical 80° Tips

SKU	Gauge
556010	23G
556011	25G



Vitreoretinal Scissors

556013

Vitreoretinal Scissors
Vertical 40° Tips

SKU	Gauge
556013	23G
556014	25G



556017

Vitreoretinal Scissors
Curved Shaft


SKU	Gauge
556017	23G
656017	25G
856017	27G



556018

Vitreoretinal Scissors
Long Curved Shaft

SKU	Gauge
556018	23G
656018	25G
856018	27G



Micro Instruments



Micro Instruments

554020

Trabeculectomy Punch

SKU	Tip Length
554020	0.75 mm
554021	1 mm
554022	1.5 mm



554045

Retinal Membrane Pick
0.5 mm Tip Length



554046

Retinal Membrane Pick
Rounded Spatula Tip,
0.5 mm Tip Length



554047


Retinal Membrane Pick
Angled 110° Tip,
0.5 mm Tip Length



Micro Instruments


554048

Retinal Membrane Pick
0.5 mm Tip Length



653059

IOL Cutting Scissors
19G Gauge



651075

IOL Cutting Forceps
23G Gauge



Micro Instruments

Micro Instruments

651072

ICL Forceps
20G Gauge



651078

ICL Forceps
Serrated, without Hole Tips, 21G Gauge



651079

ICL Forceps
Serrated, with Hole Tips, 21G Gauge




Speculums



Wire Speculums


555000

Naqahara Wire Speculum
13.5 mm Open Blades




555001

Kratz Barraquer Wire Speculum
13.5 mm Open Blades




555002

Wire Speculum
13.5 mm Open Blades



555003

Wire Speculum
13.5 mm Open Round Blades




555004

Wire Speculum
6 mm & 9 mm Open Blades



555005


Wire Speculum
9 mm Open Blades



Wire Speculums


555003

Wire Speculum
14 mm Open Blades




555007

Wire Speculum
14 mm Open Blades




555008

Wire Speculum
9 mm Open Blades



555009

Wire Speculum
11 mm Open Blades



555010

Jaffe Lid Retractors x2
16 mm Open Blades



555011

Jaffe Lid Retractors x2
12 mm Open Blades



Speculums

Wire Speculums

555012

Sprung Speculum Paediatric
6 mm Open Blades, 26.5 mm Length



555013

Sprung Speculum Neonate
4 mm Open Blades, 25 mm Length



555014

Barraquer Solid Speculum
14 mm Open Blades



555015

Barraquer Solid Speculum
14.5 mm Open Blades




555016

Wire Speculum
10 mm Open Blades



555017

Wire Speculum
14 mm Open Blades




Speculums

Wire Speculums


555018

Wire Speculum
13 mm Open Blades




555019

Wire Speculum
14 mm Open Blades



655008

Paediatric Wire Speculum
9 mm Open Blades



Adjustable Speculums

555020

Pierce Type Adjustable Speculum
16 mm Open Blades



555022

Adjustable Speculum
Swivel Open Blades



Speculums

Adjustable Speculums

555023

Lieberman Adjustable Speculum
V-Shaped Wire



555024

Lieberman Reversible Speculum
13 mm K-Wire Open Blades, 45° Angled Shaft



555025

Kershner Eye Speculum
13.5 mm Open Blades



555026

Steinert Lasik Speculum
14 mm Inner Attached Round Solid Blades



555027

Adjustable Speculum
13 mm Outer Attached Solid Square Blades



555028

Adjustable Speculum
14 mm Round Solid Blades



Adjustable Speculums

555029

Lieberman Speculum
12.5 mm Inner Attached Solid Square Blades




555030

Ginsberg Eye Speculum
13 mm Solid Blades



555031

Lindstrom-Chu Aspirating Speculum



555033

Lieberman Speculum
13 mm Open Blades



555034

Lieberman Speculum
14 mm V Type Solid Wire



655023

Adjustable Speculum
V-Wire Angled Wire



Adjustable Speculums

655035

Lieberman Speculum
Fine Open Blades



655037

Adjustable Speculum
Full Squared Blades



655038

Adjustable Speculum
Fine Full Blades



655040

Lieberman Speculum
Open Rounded Blades



655041

Adjustable Speculum
Full Rounded Blades



655042

Solid Bladed Speculum Screw Lock



Curettes



Curettes

Curettes

554030

Curette, Double Ended
Straight Blunt Tips



	Size
Tip 1	1.9 x 3.1 mm
Tip 2	2 mm

554031

Curette

SKU	Diameter
554031	1 mm
554032	1.5 mm
554033	2 mm

Retractors



Retractors

554080

Arruga Globe Retractor, Wide Retractor
10.5 mm & 2.8 mm Retractor Blades



554090

Desmarres Hook Retractor



SKU	Size
554090	Small
554091	Medium
554092	Large

555032

Lachrymal Sac Retractor
Curved Retractor Blades,
0.3 mm x 10 mm Teeth



Depressors



Depressors

Depressors

554024

Scleral Depressor
Tear Drop Shape Tip, 2 mm Size



554025

Scleral Depressor
Wide Bar, 6 mm Size



554026

Scleral Depressor
Schocket Scleral Depressor




654024

Indentor



Hooks



Hooks

554068

Jameson Strabismus Hook



SKU	Size
554068	7 mm
554069	8.5 mm

554070

Chavasse Hook

Helveston Muscle Hook, 12 mm Hook Size



554072

Retinal Detachment Hook

Oval Hold, 1.5 mm Size



554086

Gillis Hook




654070

Muscle Hook



SKU	Ending
654070	Short
654071	Long
654072	Medium

Miscellaneous



Miscellaneous

552038

Bulldog Clips
Flat Handle, 40 mm Length



552039

Blade Breaker, End Lock
Straight, Round Knurled Handle



552040

Blade Breaker, End Lock
Round Knurled Handle



MOQ

Miscellaneous

552041

Blade Handle



554049

Freer Lacrimal Chisel, Single Ends
5.8×12 mm



554074

Fine Thornton Fixation Ring
Swivel Head



SKU	Diameter	Cut Out
554074	16 mm	14 mm
554075	13 mm	11 mm

Miscellaneous

Miscellaneous

554081

Cautery, Tip, 6 mm Long Point
5×7 mm Oval Ball, Sharp Tip




651077

Foreign Body Remover



Bestsellers



Forceps

651141

Capsulorhexis Forceps
Curved Thin Shaft, Round Knurled Handle



551000

Toothed Forceps
Straight Shafts, Flat Handle, 85 mm Length,
1×2 Teeth

SKU	Size
551000	0.10 mm
551001	0.12 mm
551002	0.15 mm
551003	0.30 mm
551004	0.40 mm
551005	0.50 mm



551122

Notched Forceps
Straight Shafts, Tying Platforms,
0.25 mm Notched

SKU	Length
551122	85 mm
551123	105 mm
551124	115 mm



Forceps

551187

Colibri Toothed Forceps
Curved Shaft, 0.12 mm Size, 1x2 Teeth



551175


Tying Forceps
Curved Shaft, 6 mm Tying Platform, Flat Handle

SKU	Size
551175	85 mm
551176	105 mm
551177	115 mm



551235

Moorfields Forceps
Blunt Tipped, For Manipulating the Conjunctiva



Needle Holder

552021

Needle Holder
Curved Jaws, 9 mm Jaw Size, Round Knurled Handle

SKU	Lock
552020	Yes
552021	No



Choppers


554061

Phaco Chopper
Curved Blade, 1.8 mm Tip

Manipulators

654087

Double-Ended Manipulator
Cylinder/Y Pusher



Speculums

655035

Lieberman Speculum
Fine Open Blades



Vitreoretinal Forceps

556026

Vitreoretinal Forceps
Short End-Gripping Forceps (ILM)

SKU	Gauge
556026	23G
656026	25G
856026	27G



Bestsellers

Scissors

553014

Vannas Scissors
Curved, Sharp Pointed Tips

SKU	Tip To Pivot Length
553014	10 mm
553015	16 mm



653034

Westcott Stitch Scissors
Curved, Sharp Tips,
21 mm Tip To Pivot Length



Global Locations

Contact our HASA OPTIX distributors worldwide for premium-quality, single-use ophthalmic instruments. For an updated list of distributors, please refer to our website.



What Sets Us Apart?

Preserving the Planet with #HASARecycle

Preserving the planet for future generations is a big priority for us. This is why we are proud to have facilitated several ESG initiatives, including pioneering the development of our recycling programme #HASARecycle and creating a 10-year carbon reduction plan.

Recycling Stainless Steel Instruments

We are committed to using exclusively stainless steel in our instruments, considering its exceptional recyclability properties. With our trusted partners, we collect used instruments from hospitals and send them to recycling centres, where they are given a second life.

This recycling scheme has helped us to become shortlisted as one of Belgium's Most Sustainable Companies at the Trends Impact Awards in the Circular Economy category for two years in a row.

Decarbonising Our Supply Chain

#HASARecycle is only one of our corporate sustainability initiatives. We have also created a 10-year carbon reduction roadmap to facilitate the decarbonisation of our entire supply chain. From using mostly green energy at our manufacturing site to trialling shipments by boat instead of plane, we are dedicated to making our operations as sustainable as possible.

A standard cataract set of ours emits just 1.89kg of CO_2 - but we are committed to lowering that number to 1kg in the coming years.

B Corp Certified

Lastly, we B Corp certified, because sustainability is not just about the environment – it's about all who use it.


A B Corp certification is a designation for companies that meet high standards of social and environmental performance, accountability, and transparency. B Corps balance profit with purpose, ensuring they positively impact their employees, communities, and the environment. Through this certification, we have legally embed our commitment to social and environmental values into our business practices.


This certification gives us the opportunity to showcase our commitment to governance transparency, worker rights, and giving back to our community and customers.


Join our community to make ophthalmology more sustainable. Together, we can create a healthier future for both our patients and the planet.


Contact us for e-catalogues and discover HASA OPTIX surgical instrument solutions














info@hasaoptix.com

hasaoptix.com

+3225246388

Avenue Louise 523, 1050 Ixelles, Brussels, Belgium

Notes

Notes